

## OTA Tip-of-the-Month – November 2011

### No-Hands Retraction

William T Obrebsky MD MPH - Chief Orthopedic Trauma,  
Vanderbilt Medical Center

Sometimes a picture is worth a thousand words. This one shows how I do “Hands Free” retraction during fixation of an intertrochanteric fracture with a sliding hip screw.

A sterile towel is passed over the rod holding the vertical isolation drape. Clamped at the proper height, it suspends a towel clamp. The handle of a Bennett retractor is placed into each of the clamp’s finger loops. Tension is adjusted to position the retractors for optimal visualization, without excessive tissue pressure.

